



SeaMar® Allergy Relief

NATURAL NASAL ALLERGY PROTECTION



Category

Medical device

Form

Spray

Composition

SeaMar® - Allergy Relief solution contains:

Isotonic Sea Water solution, rich with the sea's beneficial ingredients (mineral salts and trace elements), purified water, Hyaluronic acid, and Aloe extract

Indications

SeaMar® - Allergy Relief is indicated for Allergic rhinitis (seasonal e.g. hay fever; chronic i.e. all-year-long)

Administration & Dosage

6-12 years: 1 to 2 sprays in each nostril up to 3 times daily, for 7 days or more.

>12 years: 2 to 4 sprays in each nostril up to 5 times daily, for 7 days or more.

For adults: 1 to 2 sprays in each nostril 2 to 3 times daily, for 7 days or more.

Blow nose before using **SeaMar® - Allergy Relief**. Bend neck gently forward over a washbasin and tilt your head to one side. Gently insert the nozzle into the opposite nostril, in parallel with the nasal septum, and press firmly. Turn head to the opposite side and repeat the application in the other nostril

Storage

Keep at a temperature not exceeding 30 °c in a dry place.

Keep out of reach of children

Package

Carton box contains plastic bottle with nasal applicator.





SeaMar[®] Allergy Relief

NATURAL NASAL ALLERGY PROTECTION



Benefits & Properties

- **SeaMar[®] - Allergy Relief is a 100% natural, clinically proven nasal spray that relieves symptoms associated with seasonal or chronic nasal allergies.**
- **It cleverly combines seawater isotonic solution with a unique blend of Hyaluronic acid and Aloe extract that guard the nasal tissues from contact or entry of harmful external agents (e.g. allergens, pollutants, etc), providing additional strength and protection to the nasal mucosa.**
- **SeaMar[®] - Allergy Relief relieves allergic rhinitis symptoms including blocked nose, sinus pressure, runny nose, and sneezing.**
- **SeaMar[®] - Allergy Relief cleanses the nasal cavities eliminating mucus and mucus-trapped allergens that may worsen allergic rhinitis symptoms.**